

TG-LH-EE-90-2-5

Electronics Epoxy

Description

TG-LH-EE-90-2-5 is a one-part white adhesive based on epoxy resins. It cures fast at elevated temperatures and has excellent adhesion to most pc boards and electronic components. It may be cured at 100 °C or faster at 175 °C. It has no sagging. It has a stable pot life and long shelf life even at room temperature of 2 °C. It has been designed especially for bonding electronic components. It has relatively low viscosity for faster dispense from syringes. It has high thixotropic which can control the flow of the adhesive.

Applications

Epoxy adhesive for attaching to ceramic, metals, and most plastics in electronics

Guidelines for Use

1. Thaw the epoxy to room temperature (25 °C) before use
2. Dispense the epoxy by using a syringe
3. Wipe off any excess uncured adhesive with a piece of dry cloth or tissue. Further cleaning may be achieved with tissue dabbed with iso-propanol-alcohol (IPA)
4. Cure the epoxy by heating at 100 °C for 120 minutes in a convection oven. Curing at lower temperature will require a longtime

Properties

- REACH Compliant
- ROHS Compliant

Property	TG-LH-EE-90	Unit	Test Method
Chemical type	Epoxy	-	-
Appearance	White paste	-	Visual
Mix ratio, by weight	One component	-	-
Shelf life, 2°C	6	Month	ASTM F2914
Pot life, 25°C	1	Week	ASTM F2914
Viscosity, CAP 2000+ viscometer, 25°C Cap-06@40rpm	24,000	cps	ASTM D1084
Thixotropic Index	1.8	-	-
Hardness, cured 100°C for 2 hr	90	Shore D	ASTM D2240
Shear Strength	545	kgcm ⁻²	ASTM D412
Water boil, wt gain, 100°C/1hr	0.35	%	-
Tg, DSC, cured 100°C for 2 hr	122	°C	-
Specific gravity	1.20	g/cm ³	ASTM D792
CTE, alpha-1	6.2 x 10 ⁻⁵	mm ⁻¹ °C ⁻¹	-
, alpha-2	1.8 x 10 ⁻⁴	mm ⁻¹ °C ⁻¹	-
Ionic Content, Cl	<50	ppm	-
, K	<50	ppm	-
, Na	<20	ppm	-

Packaging

5 ml EFD syringe
 10 ml EFD syringe
 30 ml EFD syringe

Recommended Cure

Temp. (°C)	Gel time	Cure time
100	12 mins	2 hrs
150	2 mins	30 mins

Storage

Store in fridge at 2°C. Tightly close original container of unused product. Storing at lower temperatures down to -20°C may prolong shelf life beyond 6 months. However it may take longer time to thaw the product

NOTICE: The information contained herein is to the best of our knowledge true and accurate. However, since the varied conditions of potential use are beyond our control, all recommendations or suggestions are presented without guarantee or responsibility on our part and users should make their own test to determine the suitability of our products in any specific situation. This product is sold without warranty either expressed or implied, of fitness for a particular purpose or otherwise, except that this product shall be of standard quality, and except to the extent otherwise stated T-Global Technology Europe and North America's invoice, quotation, or order acknowledgment. We disclaim any and liability incurred in connection with the use of information contained herein, or otherwise. All risks of such are assumed by the user. Furthermore, nothing contained herein shall be construed as a recommendation to use any process or to manufacture or to use any product in conflict with existing future patents covering any product or material or its use.